

Designing Protocols for MR imaging and Spectroscopy

Benoit Boulat
BIC – April 2005

Examples taken from the project
proposed by and worked out with

P.T. NARASIMHAN

To combine imaging and spectroscopy in
the study of the mouse brain using
intermolecular multiple quantum
coherences.

ACKNOWLEDGEMENTS

Tim Hiltner
Xiaowei Zhang

Russ Jacobs
Scott Fraser

OBJECT OF STUDY

Spectrum
Chemical signature







IMAGE → CONTRAST




TO CARE FOR OR TO DESIGN
IMAGING PROTOCOLS YOU
NEED TO KNOW THE DETAILS
OF THE TECHNIQUE THAT
GIVES RISE TO THE
CONTRAST(S).

SPECTRUM → FREQUENCIES AND AMPLITUDES



YOU NEED TO KNOW THE DETAILS OF
THE TECHNIQUE THAT GIVES RISE TO THE SPECTRUM.

SCHEMATIC OF THE BRUKER MR IMAGER



W

**WORKSTATION
(W)**

**Software
Control
And
Development**

MPS


**MAGNET
PROBES AND
SPINS (MPS)**

**Spatio-temporal
Distribution of
Magnetic fields**

Spin dynamics


MPS

MAGNETIC RESONANCE



SPIN DYNAMICS


Ultimately the dynamics imposed on the nuclear spins by the rf and the gradients will be reflected in the contrasts of the image or in the amplitudes of the resonance lines in the spectrum.



At this point the precise material nature or shape of the object of study does not count much for us. Our objective is to precisely control the nuclear spins it carries.

MPS

SPIN SPACE



- The spin manipulations are contained in the pulse sequence which is a summary of the procedures to be carried out to run the experiment. At this point the jargon used should relate to some model of spin dynamics the experimenter has in mind.

W

MPS

- Ideally we want to build a software tool that will initialize a maximum of these values with the aim of allowing the method to be used by people who may not be experts in spin dynamics.

W

- We want to optimize the method with regard to the speed of acquisition, especially if applications to living systems are sought.

W

MPS


W

What we are given or have come up with:

```

start, 10u    sgrad r2d
slice, 10u    fq8b:f1
d4    rgrad r2d
d4    groff
d3    fq1:f1
p11    ph11 ; excitation pulse
d3
d4    grad{ (t10) | (t10) | (t10) } ;coherence filter on
2m
d4    groff ;coherence filter off
d28    ;adjustment for zq
vp    ph12 ;45/135 degree pulse
d30
d4
d4    grad{ (t4) | (0) | (0) }
d11
d4    groff
;----- ECHO -----
echo, d3
d4    grad{ (t15) | (t15) | (t15) }
d5
d4    grad{ (t2) | (t2) | (t2) }
d3
p1:sp1 ph1
d3
d4    grad{ (t15) | (t15) | (t15) }
d5
d4    groff
  
```

Pulse sequence for IZQC 3D-imaging



Build a method to initialize the parameters

MPS

IZQ - FSE

IZQ
GENERATION
AND
EVOLUTION


SINGLE
QUANTUM
CONVERSION
AND
EVOLUTION

FAST SPIN
ECHO
RECORDING

Time

Generation of “unique”
contrast

To speed up
the acquisition




Fast Spin Echo (RARE)

- A method exists to initialize the parameters



MPS

FAST SPIN ECHO: SOME DETAILS



w

Fast Spin Echo



**MODULE
A
(IZQ)**

+

**MODULE
B
(FSE)**

≠

**MODULE A + MODULE B
(IZQ-FSE)**


ModuleA (IZQ) --- ModuleB (FSE)

Solutions found to solve problems of the stand alone module B do not necessarily work when some module A is placed in front of it, since the initial conditions to module B have been changed.

END of IZQ MODULE

MPS

FSE



To preserve the “unique”
contrast of IZQ we do not
want to mix M_L and M_T
inside FSE


Need to work out a
new solution for

MODULE A + MODULE B


W

W


BENCHMARK TO TEST MODIFICATIONS



Gradient along B0



Gradient orthogonal
to B0



Gradient at the magic
angle relative to B0


SPECTROSCOPY

- Spectroscopy can provide a wealth of information about living systems
- To obtain a good resolution of the spectral lines one needs to perform a good shimming of the main static field. Shimming is still pretty much an art. Only a few automated methods exist and they do not work in all situations.
- At best, shimming procedures yield a uniform field over small volumes only. Thus signal intensity is proportionately reduced.
- Field drift can be taken care of by a “field frequency lock” which is usually not available in imaging systems.
- “Block” averaging can mitigate the situation only partly and suffers other drawbacks.
- We explored the possibilities of implementing IZQ spectroscopy in our imaging system to overcome the shimming and drift problem.


MPS

Spin Systems

One spin System in a magnetic field



Two spin system in a magnetic field




Zero Quantum Coherences

Spin A – Spin X

ω_A ω_X

MPS



SQC_A: $\omega_A + \delta\omega_A(R)$

SQC_X: $\omega_X + \delta\omega_X(R)$

ZQC_{AX}: $\omega_A - \omega_X$

$B_0 + \Delta B(R)$




ZQCs are insensitive to field inhomogeneities

- One cannot directly detect through the RF coil coherences of order different from one.
- We therefore need to detect zero quantum coherence indirectly.
- We do this by using the framework of 2D spectroscopy, labeling zero quantum frequencies along the indirect dimension.

W

MPS

2D SPECTROSCOPY



W

MPS

2D IZQ SPECTROSCOPY

t_1 evolution
IZQC

Evolution under the influence
of distant dipolar field created
by the spins of the solvent

ACQ

*Insensitive to static
magnetic field
inhomogeneities*


*Generation of observable
terms*

Wide Lines


Narrow Lines

ω_1

ω_2




ONE DIMENSIONAL PROJECTIONS





IZQ spectroscopy in the live mouse brain. Dipolar evolution time = 120ms.
Relaxation delay = 5s. 4096x512 data points zero filled to 4096x1024.
Micro2.5. 20mm birdcage coil. Two averages per t1 increments.


W BENCHMARK TEST I



W BENCHMARK TEST II



W



IZQ spectroscopy in the live mouse brain. Dipolar evolution time 120ms. Coherence gradient = 20%. No water suppression. 4096x64 data points zero filled to 4096x256. Fifth order baseline correction applied in first and second dimension. Data collected 03/14/2005.